



EPP front opening isothermal container for the transport of hot, chilled and frozen meals

- Expanded polypropylene (EPP) isothermal container.
- Capacity: 83 litres.
- Lightweight and easy to handle.
- Front opening and loading container.
- Equipped with 12 runners, suitable for transporting Gastronorm 1/1 pans and submultiples.
- Perfectly stackable.
- Easy to transfer, thanks to the Maxidolly trolley.
- Designed for use with eutectic plates, to extend the transport time of hot, chilled and frozen food.
- It allows to properly operate in HACCP environments.
- The materials used for the manufacturing are suitable for food contact over the whole container surface.
- The unit's design allows accurate cleaning and sanitising (EC Reg. 852/2004 -HACCP-) including using a dishwasher.
- Fully recyclable at the end of the operating life.
- Guaranteed for use at temperatures ranging between -40° and $+120^{\circ}\text{C}$.

Model	Code	Colour	External dimensions mm	Internal dimensions mm	Capacity l	Weight Kg
Melpop AF12	MPAF0001	black	650x500xh610	540x330xh480	83	4,50
Hot GN 1/1 eutectic plate	PEGS0003	red	530x325xh30			3
Chilled GN 1/1 eutectic plate	PEGS0001	pink	530x325xh30			4
Extrachilled GN 1/1 eutectic plate	PEGS9001	sand	530x325xh30			4
Frozen GN 1/1 eutectic plate	PEGS0002	blue	530x325xh30			4
Extrafrozen GN 1/1 eutectic plate	PEGS9008	green	530x325xh30			4
Chrome handle Maxidolly trolley	MAXI0012	grey/blue mel.	525x750xh965			9,4
S.S. handle Maxidolly trolley	MAXI0011	grey/blue mel.	525x750xh965			9,4
USB Data logger	AFG002		80x33xh14			